

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SV40 Small T Antigen Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22437
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Small t antigen, ST, ST-AG

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Simian virus 40</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS