



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MMTV p10 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22491
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Gag-Pro-Pol polyprotein, Matrix protein p10, Phosphorylated protein pp21, Protein p3, Protein p8, Protein n, Capsid protein p27, Nucleocapsid protein-dUTPase, NC-dUTPase, 3.6.1.23, Protease, 3.4.23.-, Reverse transcriptase/ribonuclease H, RT, 2.7.7.49, 2.

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mouse mammary tumor virus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY