

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PTGS1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22508
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Prostaglandin G/H synthase 1, EC:1.14.99.1, Cyclooxygenase-1, COX-1, Prostaglandin H2 synthase 1, PGH synthase 1, PGHS-1, PHS 1, Prostaglandin-endoperoxide synthase 1, PTGS1, COX1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Mouse, Rat</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS