

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Streptococcus suis Enolase Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22529
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	2-phospho-D-glycerate hydro-lyase, 2-phosphoglycerate dehydratase, EC:4.2.1.11, Enolase, eno

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Streptococcus suis (strain 98HAH33)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS