

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HSD17B2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22552
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	17-beta-HSD 2, 17-beta-hydroxysteroid dehydrogenase type 2, 20 alpha-hydroxysteroid dehydrogenase, 20-alpha-HSD, E2DH, EC:1.1.1.239, EC:1.1.1.62, EDH17B2, Estradiol 17-beta-dehydrogenase 2, HSD17B2, Microsomal 17-beta-hydroxysteroid dehydrogenase, SDR9C2,

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-P37059