

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-GST 26 Nanobody (SAA1199)
PRODUCT REFERENCE	ARO-A22605
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-mFc
TARGET	Glutathione S-transferase class-mu 26 kDa isozyme, GST 26, 2.5.1.18, Sj26 antigen, SjGST

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Schistosoma japonicum</i> (Blood fluke)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Glutathione S-transferase class-mu 26 kDa isozyme, GST 26, 2.5.1.18, Sj26 antigen, SjGST, Glutathione S-transferase class-mu 26 kDa isozyme, GST 26, 2.5.1.18, Sj26 antigen, SjGST