

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human Alpha-1-antitrypsin/SERPINA1 antibody (SAA2228)
PRODUCT REFERENCE	ARO-A22608
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Alpha-1-antitrypsin/SERPINA1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Alpha-1-antitrypsin/SERPINA1, SERPINA1, Alpha-1-antitrypsin, Alpha-1 protease inhibitor, Alpha-1-antiproteinase Serpin A1, SPAAT, PRO0684, PRO2209, AAT, PI