

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human DST/Dystonin antibody (SAA2231)
PRODUCT REFERENCE	ARO-A22610
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Dystonin, 230 kDa bullous pemphigoid antigen, 230/240 kDa bullous pemphigoid antigen, Bullous pemphigoid antigen 1, BPA, Bullous pemphigoid antigen, Dystonia musculorum protein, Hemidesmosomal plaque protein, DST, BP230, BP240, BPAG1, DMH, DT, KIAA0728

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Dystonin, 230 kDa bullous pemphigoid antigen, 230/240 kDa bullous pemphigoid antigen, Bullous pemphigoid antigen 1, BPA, Bullous pemphigoid antigen, Dystonia musculorum protein, Hemidesmosomal plaque protein, DST, BP230, BP240, BPAG1, DMH, DT, KIAA07