

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human RSPO2 Recombinant Antibody (h130M23-H1L6)
PRODUCT REFERENCE	ARO-A22615
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Roof plate-specific spondin-2, R-spondin-2, hRspo2, RSPO2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Mouse</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Roof plate-specific spondin-2, R-spondin-2, hRspo2, RSPO2