

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CXCL3 Antibody (SAA2246)
PRODUCT REFERENCE	ARO-A22617
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	C-X-C motif chemokine 3 (Dendritic cell inflammatory protein 1)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-C-X-C motif chemokine 3 (Dendritic cell inflammatory protein 1), C-X-C motif chemokine 3, GRO-gamma(1-73), Growth-regulated protein gamma, GRO-gamma, Macrophage inflammatory protein 2-beta, MIP2-beta, GRO-gamma(5-73), CXCL3