

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CDKN2A/p16INK4a Nanobody (SAA2259)
PRODUCT REFERENCE	ARO-A22672
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-His
TARGET	CDKN2A/p16INK4a

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CDKN2A/p16INK4a, Tumor suppressor ARF, Alternative reading frame, ARF, Cyclin-dependent kinase inhibitor 2A, p14ARF, CDKN2A