

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Histone H1 mimotope (SSVLYGGPPSAA) Antibody (SSV mAb)
PRODUCT REFERENCE	ARO-A22695
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Histone H1 mimotope, SSVLYGGPPSAA, Histone H1.0, Histone H1', Histone H1(0), Histone H1.0 N-terminally processed, H1-0, H1F0, H1FV

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Histone H1 mimotope, SSVLYGGPPSAA, Histone H1.0, Histone H1', Histone H1(0), Histone H1.0 N-terminally processed, H1-0, H1F0, H1FV,SSVmAb