

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Cryptosporidium parvum Apical complex glycoprotein CSL Antibody (SAA2327)
PRODUCT REFERENCE	ARO-A22713
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Apical complex glycoprotein CSL, CSL (circumsporozoite-like) antigen

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Cryptosporidium parvum</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Apical complex glycoprotein CSL, CSL (circumsporozoite-like) antigen