

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CXCL4/PF4 Antibody (U.Penn. patent anti-PF4)
PRODUCT REFERENCE	ARO-A22719
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Platelet factor 4, PF-4, C-X-C motif chemokine 4, Iroplact, Oncostatin-A, Platelet factor 4 short form, Endothelial cell growth inhibitor, PF4, CXCL4, SCYB4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Platelet factor 4, PF-4, C-X-C motif chemokine 4, Iroplact, Oncostatin-A, Platelet factor 4 short form, Endothelial cell growth inhibitor, PF4, CXCL4, SCYB4