

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD138/SDC1 Antibody (587CT7.3.6.5)
PRODUCT REFERENCE	ARO-A22725
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Syndecan-1, SYND1, CD138, SDC1, SDC

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Syndecan-1, SYND1, CD138, SDC1, SDC