

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Mouse IL6R (15A7)
PRODUCT REFERENCE	ARO-A22739
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b, kappa
TARGET	Interleukin-6 receptor subunit alpha, IL-6 receptor subunit alpha, IL-6R subunit alpha, IL-6R-alpha, IL-6RA, IL-6R 1, Membrane glycoprotein 80, gp80, CD126, Soluble interleukin-6 receptor subunit alpha, sIL6R, IL6R

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Interleukin-6 receptor subunit alpha, IL-6 receptor subunit alpha, IL-6R subunit alpha, IL-6R-alpha, IL-6RA, IL-6R 1, Membrane glycoprotein 80, gp80, CD126, Soluble interleukin-6 receptor subunit alpha, sIL6R, IL6R