

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CHIKV NSP2 Antibody (Nb-A9)
PRODUCT REFERENCE	ARO-A22744
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-Fc
TARGET	Polyprotein P1234, P1234, Non-structural polyprotein, Polyprotein P123, P123, Protease nsP2, EC:3.4.22.-, EC:3.6.1.15, EC:3.6.1.74, EC:3.6.4.13, Non-structural protein 2, nsP2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Chikungunya virus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Polyprotein P1234, P1234, Non-structural polyprotein, Polyprotein P123, P123, Protease nsP2, EC:3.4.22.-, EC:3.6.1.15, EC:3.6.1.74, EC:3.6.4.13, Non-structural protein 2, nsP2