

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human PLAU/uPA Antibody (U33)
PRODUCT REFERENCE	ARO-A22765
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	Urokinase-type plasminogen activator, U-plasminogen activator, uPA, EC:3.4.21.73, Urokinase-type plasminogen activator long chain A, Urokinase-type plasminogen activator short chain A, Urokinase-type plasminogen activator chain B, PLAU

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Urokinase-type plasminogen activator, U-plasminogen activator, uPA, EC:3.4.21.73, Urokinase-type plasminogen activator long chain A, Urokinase-type plasminogen activator short chain A, Urokinase-type plasminogen activator chain B, PLAU