

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HTLV-1 Glycoprotein 46/gp46 Antibody (SAA2372)
PRODUCT REFERENCE	ARO-A22778
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Envelope glycoprotein gp62, Env polyprotein, Surface protein, SU, Glycoprotein 46, gp46, Transmembrane protein, TM, Glycoprotein 21, gp21, env

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human T-cell leukemia virus 1 (HTLV-1)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Envelope glycoprotein gp62, Env polyprotein, Surface protein, SU, Glycoprotein 46, gp46, Transmembrane protein, TM, Glycoprotein 21, gp21, env