

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse CD11c/ITGAX Antibody (N418)
PRODUCT REFERENCE	ARO-A22787
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Integrin alpha-X, CD11 antigen-like family member C, Leukocyte adhesion glycoprotein p15095 alpha chain, Leukocyte adhesion receptor p15095, CD11c, Itgax

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Integrin alpha-X, CD11 antigen-like family member C, Leukocyte adhesion glycoprotein p15095 alpha chain, Leukocyte adhesion receptor p15095, CD11c, Itgax