

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human VMAT2 Antibody (JP1217)
PRODUCT REFERENCE	ARO-A22803
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Synaptic vesicular amine transporter, Solute carrier family 18 member 2, Vesicular amine transporter 2, VAT2, Vesicular monoamine transporter 2, SLC18A2, SVMT, VMAT2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Pig</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Synaptic vesicular amine transporter, Solute carrier family 18 member 2, Vesicular amine transporter 2, VAT2, Vesicular monoamine transporter 2, SLC18A2, SVMT, VMAT2