

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-2,4-Dichlorophenoxyacetic acid/2,4-D Antibody (SAA2414)
PRODUCT REFERENCE	ARO-A22852
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	2,4-dichlorophenoxyacetic acid, CAS: 94-75-7, 2-(2,4-dichlorophenoxy)acetic acid, 2,4-D, (2,4-Dichlorophenoxy)acetic acid

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-2,4-dichlorophenoxyacetic acid, CAS: 94-75-7, 2-(2,4-dichlorophenoxy)acetic acid, 2,4-D, (2,4-Dichlorophenoxy)acetic acid