

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human Heparanase Neutralizing Antibody (A54)
PRODUCT REFERENCE	ARO-A22865
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	Heparanase, EC:3.2.1.166, Endo-glycosaminidase, Heparanase-1, Hpa1, Heparanase 8 kDa subunit, Heparanase 50 kDa subunit, HPSE, HEP, HPA, HPA1, HPR1, HPSE1, HSE1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Heparanase, EC:3.2.1.166, Endo-glycosaminidase, Heparanase-1, Hpa1, Heparanase 8 kDa subunit, Heparanase 50 kDa subunit, HPSE, HEP, HPA, HPA1, HPR1, HPSE1, HSE1