

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human S100A4 Antibody (5C3)
PRODUCT REFERENCE	ARO-A22871
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Protein S100-A4, Calvasculin, Metastasin, Placental calcium-binding protein, Protein Mts1, S100 calcium-binding protein A4, S100A4, CAPL, MTS1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Mouse, Cynomolgus monkey, Rhesus macaque</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Protein S100-A4, Calvasculin, Metastasin, Placental calcium-binding protein, Protein Mts1, S100 calcium-binding protein A4, S100A4, CAPL, MTS1