



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT NAME</b>      | Anti-Human CD328/SIGLEC7 Antibody (SAA2438)  |
| <b>PRODUCT REFERENCE</b> | ARO-A22885   |
| <b>CLONALITY</b>         | Monoclonal Antibody  |
| <b>ISOTYPE</b>           | IgG1, kappa  |
| <b>TARGET</b>            | Sialic acid-binding Ig-like lectin 7, Siglec-7, Adhesion inhibitory receptor molecule 1, AIRM-1, CDw328, D-siglec, QA79 membrane protein, p75, CD328, SIGLEC7, AIRM1 |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>ORIGIN SPECIES</b>          | <i>Homo sapiens (Human)</i>  |
| <b>FORM</b>                    | Liquid   |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4)   |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-Sialic acid-binding Ig-like lectin 7, Siglec-7, Adhesion inhibitory receptor molecule 1, AIRM-1, CDw328, D-siglec, QA79 membrane protein, p75, CD328, SIGLEC7, AIRM1  |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**