

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human F9/Factor IX Antibody (SAA2440)
PRODUCT REFERENCE	ARO-A22887
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4, kappa
TARGET	Coagulation factor IX, EC:3.4.21.22, Christmas factor, Plasma thromboplastin component, PTC, Coagulation factor IXa light chain, Coagulation factor IXa heavy chain, F9

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Coagulation factor IX, EC:3.4.21.22, Christmas factor, Plasma thromboplastin component, PTC, Coagulation factor IXa light chain, Coagulation factor IXa heavy chain, F9