

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human IA-2/PTPRN Antibody (96/3)
PRODUCT REFERENCE	ARO-A22935
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	ICA 512, ICA3, ICA512, ICA512-CCF, ICA512-N-terminal fragment, ICA512-NTF, ICA512-TMF, ICA512-cleaved cytosolic fragment, ICA512-transmembrane fragment, Islet cell antigen 512, Islet cell autoantigen 3, PTP IA-2, PTPRN, R-PTP-N, Receptor-type tyrosine-pro

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Mouse, Rat</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ICA 512, ICA3, ICA512, ICA512-CCF, ICA512-N-terminal fragment, ICA512-NTF, ICA512-TMF, ICA512-cleaved cytosolic fragment, ICA512-transmembrane fragment, Islet cell antigen 512, Islet cell autoantigen 3, PTP IA-2, PTPRN, R-PTP-N, Receptor-type tyrosin