

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CLEC6A/Dectin-2 Antibody (BDC-3042)
PRODUCT REFERENCE	ARO-A22938
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	C-type lectin domain family 6 member A, C-type lectin superfamily member 10, CLEC6A, CLECSF10, DC-associated C-type lectin 2, DECTIN2, Dectin-2, Dendritic cell-associated C-type lectin 2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-C-type lectin domain family 6 member A, C-type lectin superfamily member 10, CLEC6A, CLECSF10, DC-associated C-type lectin 2, DECTIN2, Dectin-2, Dendritic cell-associated C-type lectin 2