

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-GPA33 & CD3E Bispecific Antibody (MGD007)
PRODUCT REFERENCE	ARO-A22939
CLONALITY	Monoclonal Antibody
ISOTYPE	Dart-Fc
TARGET	Cell surface A33 antigen, GPA33, Glycoprotein A33, CD3E, CD3e, T-cell surface antigen T3/Leu-4 epsilon chain, T-cell surface glycoprotein CD3 epsilon chain, T3E

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Cell surface A33 antigen, GPA33, Glycoprotein A33, CD3E, CD3e, T-cell surface antigen T3/Leu-4 epsilon chain, T-cell surface glycoprotein CD3 epsilon chain, T3E