

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-AAV9 VP1/Capsid protein VP1 Antibody (SAA2503)
PRODUCT REFERENCE	ARO-A22945
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	capsid protein VP1 [Adeno-associated virus], cap, clade F, Phospholipase A2-like domain; pfam08398, Parvovirus coat protein VP2; pfam00740

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Adeno-associated virus 8 (AAV8)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-capsid protein VP1 [Adeno-associated virus], cap, clade F, Phospholipase A2-like domain; pfam08398, Parvovirus coat protein VP2; pfam00740