

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD182/CXCR2 Antibody (SAA1360)
PRODUCT REFERENCE	ARO-A22956
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-Fc
TARGET	C-X-C chemokine receptor type 2, CD182, CDw128b, CXC-R2, CXCR-2, CXCR2, GRO/MGSA receptor, High affinity interleukin-8 receptor B, IL-8 receptor type 2, IL-8R B, IL8RB

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-C-X-C chemokine receptor type 2, CD182, CDw128b, CXC-R2, CXCR-2, CXCR2, GRO/MGSA receptor, High affinity interleukin-8 receptor B, IL-8 receptor type 2, IL-8R B, IL8RB