

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human IL38/IL1F10 Antibody (SAA2514)
<b>PRODUCT REFERENCE</b>	ARO-A22960
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	FIL1 theta, FIL1T, Family of interleukin 1-theta, IL-1 theta, IL-1F10, IL-1HY2, IL-38, IL1F10, IL1HY2, IL38, Interleukin-1 HY2, Interleukin-1 family member 10, Interleukin-1 theta, Interleukin-38

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Human, Mouse</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-FIL1 theta, FIL1T, Family of interleukin 1-theta, IL-1 theta, IL-1F10, IL-1HY2, IL-38, IL1F10, IL1HY2, IL38, Interleukin-1 HY2, Interleukin-1 family member 10, Interleukin-1 theta, Interleukin-38