

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse BCMA & CD3E Bispecific Antibody (SAA2527)
PRODUCT REFERENCE	ARO-A22976
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	RecName: Full=Tumor necrosis factor receptor superfamily member 17; AltName: Full=B-cell maturation protein; AltName: CD_antigen=CD269, Tnfrsf17, B-cell maturation protein, Tumor necrosis factor receptor superfamily member 17. /id=PRO_0000058936., Tumor n

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-RecName: Full=Tumor necrosis factor receptor superfamily member 17; AltName: Full=B-cell maturation protein; AltName: CD_antigen=CD269, Tnfrsf17, B-cell maturation protein, Tumor necrosis factor receptor superfamily member 17. /id=PRO_0000058936., Tu