

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human EGFR/ERBB1/HER1 Antibody (SAA2536)
PRODUCT REFERENCE	ARO-A22987
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv-His
TARGET	EC:2.7.10.1, EGFR, ERBB, ERBB1, Epidermal growth factor receptor, HER1, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-EC:2.7.10.1, EGFR, ERBB, ERBB1, Epidermal growth factor receptor, HER1, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1