

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Bet v 1/Allergen Bet v I-A Nanobody (SAA0668)
PRODUCT REFERENCE	ARO-A22997
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-His
TARGET	Allergen Bet v I-A, BETVI, BETVIA, Bet v 1-A, Major pollen allergen Bet v 1-A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Betula pendula (European white birch) (Betula verrucosa)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Allergen Bet v I-A, BETVI, BETVIA, Bet v 1-A, Major pollen allergen Bet v 1-A