

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-HPV16 L1/Major capsid protein L1 Antibody (SAA2553)
PRODUCT REFERENCE	ARO-A22999
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	L1, Major capsid protein L1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human papillomavirus type 16 (HPV16)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-L1, Major capsid protein L1