



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-E3G Antibody (SAA2650) |
| PRODUCT REFERENCE | ARO-A23079 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2, kappa |
| TARGET | EI3G, Estriol 3-O-(beta-D-glucuronide), Estriol 3-(beta-D-glucuronide), Estriol-3-glucuronide E3G, Estrone beta-D-glucuronide, Estrone 3-O-(beta-D-glucuronide), Estrone-3-glucuronide |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| ORIGIN SPECIES | <i>Mus musculus (Mouse)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-EI3G, Estriol 3-O-(beta-D-glucuronide), Estriol 3-(beta-D-glucuronide), Estriol-3-glucuronide E3G, Estrone beta-D-glucuronide, Estrone 3-O-(beta-D-glucuronide), Estrone-3-glucuronide |

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY