

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human Human CYP19A1 Antibody (SAA2693)
PRODUCT REFERENCE	ARO-A23092
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	ARO1, Aromatase, CYAR, CYP19, CYP19A1, CYPXIX, Cytochrome P-450AROM, Cytochrome P450 19A1, EC:1.14.14.14, Estrogen synthase

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ARO1, Aromatase, CYAR, CYP19, CYP19A1, CYPXIX, Cytochrome P-450AROM, Cytochrome P450 19A1, EC:1.14.14.14, Estrogen synthase