

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-SNAP-25 (BoTox-A Cleaved) Antibody (SAA2659)
PRODUCT REFERENCE	ARO-A23095
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b
TARGET	SNAP, SNAP-25, SNAP25, SUP, Super protein, Synaptosomal-associated 25 kDa protein, Synaptosomal-associated protein 25

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-SNAP, SNAP-25, SNAP25, SUP, Super protein, Synaptosomal-associated 25 kDa protein, Synaptosomal-associated protein 25