

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse CCL20/MIP-3-alpha Antibody (SAA2679)
PRODUCT REFERENCE	ARO-A23116
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Beta-chemokine exodus-1, C-C motif chemokine 20, CC chemokine LARC, CC chemokine ST38, Ccl20, Larc, Liver and activation-regulated chemokine, MIP-3-alpha, Macrophage inflammatory protein 3 alpha, Scya20, Small-inducible cytokine A20

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Beta-chemokine exodus-1, C-C motif chemokine 20, CC chemokine LARC, CC chemokine ST38, Ccl20, Larc, Liver and activation-regulated chemokine, MIP-3-alpha, Macrophage inflammatory protein 3 alpha, Scya20, Small-inducible cytokine A20