

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Hepatitis E virus/HEVHEV Capsid protein/ORF2 Antibody (SAA2698)
PRODUCT REFERENCE	ARO-A23132
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	ORF2s, Pro-secreted protein ORF2, Protein ORF2, Secreted protein ORF2, pORF2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Hepatitis E virus genotype 1 (isolate Human/Pakistan/Sar-55) (HEV-1)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ORF2s, Pro-secreted protein ORF2, Protein ORF2, Secreted protein ORF2, pORF2