



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human TGM2 Antibody (SAA2703)
<b>PRODUCT REFERENCE</b>	ARO-A23135
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgGA, kappa
<b>TARGET</b>	EC:2.3.1.-, EC:2.3.2.13, EC:3.4.-.-, EC:3.5.1.44, Erythrocyte transglutaminase, G(h), Heart G alpha(h), Isopeptidase TGM2, Protein G alpha(h), Protein-glutamine deamidase TGM2, Protein-glutamine dopaminytransferase TGM2, Protein-glutamine gamma-glutamylt

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-EC:2.3.1.-, EC:2.3.2.13, EC:3.4.-.-, EC:3.5.1.44, Erythrocyte transglutaminase, G(h), Heart G alpha(h), Isopeptidase TGM2, Protein G alpha(h), Protein-glutamine deamidase TGM2, Protein-glutamine dopaminytransferase TGM2, Protein-glutamine gamma-glut

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**