

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human Frataxin Antibody (SAA2739)
<b>PRODUCT REFERENCE</b>	ARO-A23164
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-Fc
<b>TARGET</b>	EC:1.16.3.1, Extramitochondrial frataxin, FRDA, FXN, Frataxin intermediate form, Frataxin mature form, Frataxin(56-210), Frataxin(78-210), Frataxin(81-210), Frataxin, mitochondrial, Friedreich ataxia protein, Fxn, X25, d-FXN, i-FXN, m56-FXN, m78-FXN, m81-

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-EC:1.16.3.1, Extramitochondrial frataxin, FRDA, FXN, Frataxin intermediate form, Frataxin mature form, Frataxin(56-210), Frataxin(78-210), Frataxin(81-210), Frataxin, mitochondrial, Friedreich ataxia protein, Fxn, X25, d-FXN, i-FXN, m56-FXN, m78-FXN,