

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD261/TNFRSF10A/DR4 Antibody (SAA2670)
PRODUCT REFERENCE	ARO-A23213
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	APO2, CD261, DR4, Death receptor 4, TNF-related apoptosis-inducing ligand receptor 1, TNFRSF10A, TRAIL receptor 1, TRAIL-R1, TRAILR1, Tumor necrosis factor receptor superfamily member 10A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-APO2, CD261, DR4, Death receptor 4, TNF-related apoptosis-inducing ligand receptor 1, TNFRSF10A, TRAIL receptor 1, TRAIL-R1, TRAILR1, Tumor necrosis factor receptor superfamily member 10A