

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ATP5PB Monoclonal Antibody (1A480)
PRODUCT REFERENCE	ARO-A23251
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	ATP synthase peripheral stalk subunit b mitochondrial, ATP synthase F(0) complex subunit B1 mitochondrial, ATP synthase peripheral stalk-membrane subunit b, ATP synthase proton-transporting mitochondrial F(0) complex subunit B1, ATP synthase subunit b, AT

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ATP synthase peripheral stalk subunit b mitochondrial, ATP synthase F(0) complex subunit B1 mitochondrial, ATP synthase peripheral stalk-membrane subunit b, ATP synthase proton-transporting mitochondrial F(0) complex subunit B1, ATP synthase subunit