

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human INPP5D Monoclonal Antibody (1A540)
PRODUCT REFERENCE	ARO-A23308
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	Phosphatidylinositol 345-trisphosphate 5-phosphatase 1, EC:3.1.3.86, Inositol polyphosphate-5-phosphatase D, EC:3.1.3.56, Inositol polyphosphate-5-phosphatase of 145 kDa, SIP-145, Phosphatidylinositol 45-bisphosphate 5-phosphatase, EC:3.1.3.36, SH2 domain

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Phosphatidylinositol 345-trisphosphate 5-phosphatase 1, EC:3.1.3.86, Inositol polyphosphate-5-phosphatase D, EC:3.1.3.56, Inositol polyphosphate-5-phosphatase of 145 kDa, SIP-145, Phosphatidylinositol 45-bisphosphate 5-phosphatase, EC:3.1.3.36, SH2 d