

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human F8 Monoclonal Antibody (1A603)
PRODUCT REFERENCE	ARO-A23368
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	Coagulation factor VIII, Antihemophilic factor, AHF, Procoagulant component, Factor VIIIa heavy chain 200 kDa isoform, Factor VIIIa heavy chain 92 kDa isoform, Factor VIII B chain, Factor VIIIa light chain, F8

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Coagulation factor VIII, Antihemophilic factor, AHF, Procoagulant component, Factor VIIIa heavy chain 200 kDa isoform, Factor VIIIa heavy chain 92 kDa isoform, Factor VIII B chain, Factor VIIIa light chain, F8