

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Staphylococcus aureus srtA/Sortase A Protein, N-His-SUMO
PRODUCT REFERENCE	ARO-P10255
ORIGIN SPECIES	Staphylococcus aureus (strain NCTC 8325 / PS 47)
SOURCE	
FRAGMENT	Gln60-Lys206
TAG	N-His-SUMO Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	28.84 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.