

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Recombinant CCHFV GP/Envelope glycoprotein Protein, N-His
<b>PRODUCT REFERENCE</b>	ARO-P10283
<b>ORIGIN SPECIES</b>	Crimean-Congo hemorrhagic fever virus (strain Nigeria/IbAr10200/1970) (CCHFV)
<b>SOURCE</b>	
<b>FRAGMENT</b>	Ser520-Ser698
<b>TAG</b>	N-His Tag

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>PRODUCTION HOST</b>	<i>Escherichia coli (E. coli)</i>
<b>MOLECULAR WEIGHT</b>	21.96 kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.