

## PRODUCT DATASHEET

### OVERVIEW

---

**PRODUCT NAME** Recombinant Rat VEGFA Protein, N-His

**PRODUCT REFERENCE** ARO-P10887

**ORIGIN SPECIES** Rattus norvegicus (Rat)

**SOURCE**

**FRAGMENT** Ala27-Arg214

**TAG** N-His Tag

### PRODUCT INFORMATION

---

**EXPRESSION SYSTEM** Prokaryotic expression

**PRODUCTION HOST** *Escherichia coli* (*E. coli*)

**MOLECULAR WEIGHT** 24.37 kDa

**FORM** Lyophilized

**BUFFER** Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

**STABILITY & STORAGE** Store at 4°C for short term while lyophilized.  
Freeze and store at -20°C or -80°C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Avoid freeze/thaw cycles, aliquot the product according to your needs.